

Collecting Grain Samples:

1. You will require plastic sample bags and a pen to label bags. Wash your hands thoroughly before sampling or wear a clean rubber glove.
2. Collect a minimum of 10 sub-samples (approximately 1 hand full of grain) and mix. Mark identification details. Fill in this sheet as your identification details will be the ones we use when returning your results.
3. Store samples in a cool place such as an esky or refrigerator between sampling and posting. Post/freight samples ASAP ensuring all details are completed. We require at least 10 working days to carry out germinations.

| | | | |
|------------------------|--|----------------------|--|
| BILLING NAME: | | | |
| POSTAL ADDRESS: | | | |
| CONTACT NAME: | | | |
| PHONE: | | FAX: | |
| email: | | | |
| PROPERTY NAME: | | DATE SAMPLED: | |
| Any Comments: | | | |

| | 1 | 2 |
|---|----------|----------|
| Seed Type: | | |
| Sample Variety: | | |
| Id/Line number: | | |
| Test Requested: <i>(eg Germination, purity, protein):</i> | | |
| Special Conditions required? | | |
| Comments: | | |